

***aquapick***<sup>®</sup>

www.aquapick.com

| **Aquapick Co., Ltd.** |

20F, SK V1 Center, 171, Gasan digital 1-ro, Geumcheon-gu, Seoul, Korea | **E-mail.** info@aquapick.com

Customer Service: **+82-1588-2809**

***aquapick***<sup>®</sup>

Officially recommended products by  
Korean Dental Association

2021 AQUAPICK CATALOGUE

*“ A leader  
in oral care solutions „*








# aquapick®


*Company Profile*



Lobby



Book Cafe



Office Area



Office Area




Product Demonstration Room



Product Display Room



VIP Reception Room




Conference Hall




Office Area



Island Kitchen



Cafeteria



Aquapick Award and Certification Display





## *Message from CEO*

Although improved living standards have enhanced awareness about the care of the teeth and gums, proper care is still not easy to implement in everyday life. This is primarily due to the presence of blind spots where brushing alone does not offer a good solution. Even brushing three times a day still fails to completely prevent small and large problems with the gums and teeth effectively. These blind spots are areas where a brush cannot reach. Pieces of food left around these areas become hotbeds for growing Plaque and accelerating the formation of tartar, presenting a threat to the health of the teeth and gums.

Aquapick, which specializes in developing and manufacturing oral care systems, has made efforts to protect and improve your oral health by offering simple and easy solutions to oral care, which requires more than brushing. Thanks to these efforts, Aquapick has won the Presidential Citation which was given by the President of Korea. This is given to the best of best companies that have contributed to social and economic development. Aquapick has taken the lead in the oral care market and promises to continue to afford customers greater satisfaction with better quality and service through intensive investment in R&D and continuous product innovation.

Aquapick will not be satisfied with today's status quo as the best oral care system in Korea but will continue to make every effort to fulfill its social responsibility by expanding the base of the oral care system and developing diverse products, enabling more people to take care of their oral health with ease.

Furthermore, Aquapick will pursue diverse activities for the enhancement of customer satisfaction including fulfillment of thorough warranty services. Aquapick has grown with your continued encouragement. We are thankful for your continued faith and trust and we will press forward in our efforts to meet and exceed your expectations with the highest quality products and services.

Thank you.

CEO of Aquapick Co., Ltd.

A handwritten signature in black ink that reads "Kevin Lee".







# *aquapick*<sup>®</sup>

The products that we are selling is not at their highest quality. It is because we are making better products at the moment.

We are not the 'yes man.'

We are the 'no man.'

We always say 'no.'

We will always say "no" until you say "yes."

Fulfilling customer satisfaction, we will pay back our gratitude with greater customer service.

\* 1 year warranty





Must-have for  
your oral health  
**Aquapick!**



# GOOD HABIT

Correct habits for strong teeth and healthy gums



Removes food debris while  
massaging gums gently



Prevents peri-implantitis



Ideal for orthodontic patients



Takes care of pregnant  
women's teeth (before & after)



Prevention of  
kids' tooth decay



Remove bad breath


**Aquapick** makes teeth and gums healthy and strong,  
by cleaning with a pulsating stream of water at 2,200 pulses/minute.





# HEALTH

Oral health is the beginning of general health.



→ If you clean your teeth carelessly,



Poor oral hygiene accelerates aging process and reduces cognitive ability.

- Journal of American gerontology (Denmark Research Team) -



Greater risks for pneumonia, asthma, and laryngopharyngitis.

-February, 2002, Issue of the British Journal on the Periodontal Membrane-



Greater risk (up to 2.8 times greater than healthy people) for cardiac arrest.

-International Dental Research Association -



Greater risk (up to 7.8 times greater than healthy people) for premature birth.

-International Dental Research Association -



Infection or inflammation in the roots of teeth, may increase the risk of diabetes by lowering insulin levels and cause difficulty in controlling blood sugar levels.

- Taiwanese Dental Health Association-

# RESEARCH


Aquapick is an oral health care device with clinically verified results.

## Comparison of interdental plaque using Aquapick

Clinical Institution: The Korean Academy of Oral Health  
(done by Paired-test)  
Researcher: Dr. Bae  
(School of dentistry, Seoul National University)

Aquapick is an oral health care device with clinically verified results.

Interdental plaque index was drastically reduced by using Aquapick.




## Comparison between Aquapick, dental floss and interdental brush.

Clinical Institution: The Korean Academy of Oral Health  
(done by one-way analysis of variance)  
Standard deviation 0.67~1.21  
Researcher: Dr. Bae  
(School of dentistry, Seoul National University)

Aquapick is an oral health care device with high functionality, convenience and satisfaction.

Aquapick has outstanding functionality and satisfaction after use. It is also very convenient when you clean your teeth that are difficult for dental floss or interdental brush to reach



Aquapick is helpful to prevent gum disease.

**aquapick**







# *aquapick*®

Aquapick is a daily **must-have**  
**device** for modern people.



| AQ-230



| AQ-120



| AQ-350

**NEW**


# Aquapick Oral Irrigator AQ-350



Officially recommended  
by the Korean Dental Association

Product Name: Aquapick Oral Irrigator  
Model: AQ-350  
Color: Pink/Gray  
Components: Main Body, Water Tank, Jet Tip 4EA,  
Jet Tip Holder 4EA  
Electrical Currency: 100-230V~, 50/60Hz  
Power Consumption: 24W  
Dimension: 150x90x210mm (WxDxH)  
Water Tank Capacity: 600ml  
Pulsations: 2,400 pulses/minute

Jet Tip Color: ● ● ● ●



| Pulsating stream of water  
at 2,400 pulses/minute

Effective on removal of food  
between tooth & gums




| Convenient handle  
mounting with  
ergonomic design

Effective on removal of food  
between tooth & gums



| AUTO START

If you pull up the handle, it  
**automatically starts**  
**functioning** by activating  
the sensor.  
(The world's first development)



| 600ml water tank

Sufficiently usable volume





# Aquapick Oral Irrigator AQ-310



Product Name: Aquapick Oral Irrigator  
Model: AQ-310  
Components: Main Body, Water Tank,  
Jet Tip 4EA, Special tips 4EA  
(Tongue scrapper,  
Orthodontic tip, Pocket tip)  
Electrical Currency: 100-230V~, 50/60Hz  
Power Consumption : 18W  
Dimension: 140x115x200mm (WxDxH)  
Water Tack Capacity: 600ml  
Pulsations: 1,600 pulses/min

Jet Tip Color: ● ● ● ●

## Sophisticated | water pressure controller

Men and women of all ages  
can use it, thanks to  
10 steps of sophisticated  
water pressure controller.



## Jet tip case |

Easily store jet tips  
with jet tip case.



## 600ml water tank |

Sufficiently usable volume



## | 3 special tips

Three special tips for  
perfect oral care allow  
three-dimensional mouthwash



## | 4 different colors of clear jet tip

All family members can use it



# Aquapick Oral Irrigator AQ-300



Product Name: Aquapick Oral Irrigator  
Model: AQ-300  
Components: Main Body, Water Tank, Jet Tip 4EA,  
Jet Tip Holder 4EA  
Electrical Currency: 100-230V~, 50/60Hz  
Power Consumption: 24W  
Dimension: 150x90x210mm (WxDxH)  
Water Tank Capacity: 600ml  
Pulsations: 2,200 pulses/minute

Jet Tip Color: ● ● ● ●



## | Pulsating stream of water at 2,200 pulses/minute

Effective on removal of food between  
teethes&gums and massage on the gums



## | Water drop shaped tip holders

Sanitarily store jet tips  
(can be used as tooth brush holders)



## | 4 different colors of clear jet tip

All family members can use it



## | Easy control of water pressure

Men and women of all ages can use it  
due to sophisticated water pressure control





**NEW**

# Aquapick Portable Oral Irrigator AQ-240




Product Name: Aquapick Portable  
Oral Irrigator  
Model: AQ-240  
Components: Main Body, Water Tank,  
Jet Tip 4EA, Charger Cable  
Electrical Current: 5V  
Power Consumption: 3W  
Dimension: 68x48x132mm (WxDxH)  
Weight: 422g  
Water Tank capacity: 450ml  
Pulsations: 1,400 pulses/min  
2 cleaning modes: Clean, Massage

Jet Tip Color: ● ● ● ●



• **Pouch for  
Traveling**

Pouch for storing and  
protecting the product  
when traveling



• **Large Capacity  
Water Tank**

Sufficiently use the product  
with 450ml water  
tank capacity



• **Water Pressure  
Control**

Adjust the water pressure in 5 levels  
by moving the water pressure control  
button up and down



• **Portable Size**

Neat and small size,  
suitable for carrying around  
when traveling



# Aquapick Cordless Oral Irrigator AQ-230



Officially recommended  
by the Korean Dental Association

Product Name: Aquapick Cordless Oral Irrigator  
Model: AQ-230  
Components: Handle, Water Tank, Charger  
Jet Tip 2EA  
Electrical Currency: 100-230V~, 50/60Hz  
Power Consumption: 24W  
Dimension: 60x85x275mm (WxDxH)  
Weight: 280g  
Water Tank capacity: 200ml  
Pulsations: 1,400 pulses/min  
3 cleaning mode: Clean, Soft, Massage

Jet Tip Color: ● ●



- Bigger sized water tank  
Expanded capacity of **200ml**  
(Normal water tank: 150ml)



- 3 cleaning mode  
**Clean/Soft/Massage**  
Available for  
customized management



- Conveniently Portable  
Optimal size for  
**carrying around**

\* Case separately  
purchasable






# Aquapick Cordless Oral Irrigator AQ-220



Product Name: Aquapick Cordless Oral Irrigator  
Model: AC-220  
Components: Handle, Water Tank, Charger  
Jet Tip 2EA  
Electrical Currency: 110-240V~, 50/60Hz  
Power Consumption: 0.7-1.7W  
Dimension: 60x75x275mm (WxDxH)a  
Weight: 260g  
Water Tank Capacity: 150ml  
Pulsations: 1,400 pulses/minute  
3 cleaning mode: Normal, Soft, Pulse

Jet Tip Color: ● ●



- 3 Cleaning mode enables customized management
  - Normal** : General cleaning mode
  - Soft** : Soft cleaning mode for sensitive teeth
  - Pulse** : Gum massaging mode
  - CH** : Charging Indicator



- Upgraded charging method  
**non-contact charging method!**  
It will start charging just by placing on the spot.  
Easy & safe!



Takes 24 hours  
to fully charge



The battery can last up to 2 weeks  
when fully charged

NEW

# Aquapick Cordless Oral Irrigator AQ-205



Product Name: Aquapick Cordless Oral Irrigator

Model: AQ-205

Components: Handle, Water Tank,  
Charging Cable, Jet Tip 2EA

Electrical Currency: 3.7V

Power Consumption: 5W

Dimension: 58x40.8x211mm (WxDxH)

Weight: 369g

Water Tank capacity: 170ml

Pulsations: 1,700 pulses/minute

3 cleaning modes: Clean, Soft, Massage

Jet Tip Color: ● ●



## • Easy to Carry

There is a **jet tip hook** inside the water tank at the bottom of the main body, so it can be easily carried.



## • Convenient magnetic charging

USB type cable and **magnetic charging method** make it more **convenient** and **safer** to use than general charging cables.



## • Large capacity battery

It can be used for up to **30 days** without charging when it is fully charged for **4 hours**.







# Aquapick Sonic Electric Toothbrush AQ-120




Product Name: Aquapick Sonic Electric Toothbrush  
Model: AQ-120  
Components: Handle, Charger,  
Brush Head 2EA  
Electrical Currency: 100~230V / 50~60Hz  
Power consumption: 24W  
Dimension: 30x30x250mm (WxDxH)  
Weight: 120g  
Vibration: 40,000 pulses/min  
3 Brushing Mode: Clean, Soft, Massage


Brush Head Color: ● ●



- **Portable size**  
Convenient to carry around because of the **small size and light weight**



- **Smart Interval Timer**  
Smart tooth brush which pauses in **every 30 seconds** It recommends proper brushing habit for children



- **3 different brushing mode**  
Clean : General mode for brushing teeth  
Soft : Washes sensitive spots or teeth  
Massage: Massages gums softly



- **Non-contact charging method**  
Designed to be **safe with waterproof** level of IPX 7



\* Case separately purchasable

NEW

# Aquapick Sonic Electric Toothbrush AQ-102



Product Name: Aquapick Sonic  
Electric Toothbrush

Model: AQ-102

Color: Midnight Black/Noble white/  
Pink Blossom/Mystic Red

Componets: Handle, Charger, Brush Head 2EA

Electrical Currency: 3.7V

Power Consumption : 2W

Vibration: 36,000 pulses/min

Dimension: 29.5x29.5x245mm (WxDxH)

Weight: 105g


2 Brushing Modes: Clean, Sensitive

Brush Head Color: ● ●



- Sonic vibration bubble with more than 30,000 times/ min

Gets rid of food debris which is hidden here and there by using sonic waves.



- Smart Interval Timer

Smart tooth brush which pauses in every 30 seconds It recommends proper brushing habit for children



- USB charging type

Convenient to carry around due to small size and light weight





NEW

# Aquapick Sonic Electric Toothbrush AQ-101



Product Name: Aquapick Sonic  
Electric Toothbrush

Model: AQ-101

Componets: Handle, Charger, Brush Head 2EA

Electrical Currency: 3.7V

Power Consumption : 1.85W

Vibration: 40,000 pulses/min

Size: 21x21x224mm (WxDxH)

Weight: 88g


3 Brushing Modes: Clean, Soft, Massage

Brush Head Color: ● ●



- Dupont company's (U.S.) brush

Minimize damage on teeth enamel by using Dupont company's brush which has the highest quality in the world.



- Smart Interval Timer

Smart tooth brush which pauses in **every 30 seconds** It recommends proper brushing habit for children



- USB charging type

**Small sized and light weighted**



# Aquapick Sonic Pulse Electric Toothbrush AQ-110



Product Name: Aquapick Sonic Pulse  
Electric Toothbrush

Model: AQ-110

Components: Handle, Charger,  
Brush Head 3EA

Electrical Currency: 220V~ / 60Hz

Power Consumption: 2.5W

Dimension: 30x30x250mm (WxDxH)

Weight: 200g


Vibration: 13,000~15,000 times/minute

3 Brushing Mode: Normal, Soft, Pulse

Brush Head Color: ○ ● ●



- **Speedy Sonic Vibration**  
Approximately **13,000~15,000 times per minute**



- **Interval Timer of 30 seconds**  
Smart tooth brush which stops **every 30 seconds**



- **3 Brush heads with ergonomic design**


**a Professional Brush Head**  
2 brush bristles are designed to be integrated and to wash tooth surface and space between teeth and gums at once.

**b Slimsonic Brush Head**  
Even weak teeth and gums can be gently brushed with ergonomically designed brush head.

**c Inter-dental Brush Head**  
Suitable for washing orthodontic devices, crowns and teeth of malocclusion

- **3 different brushing mode**

| **Normal** : General mode for brushing teeth .....  
| **Soft** : Washes sensitive spots .....  
| **Pulse** : Massages gums softly .....



- **Special UV sterilizer**

99.9% bactericidal effect against Escherichia coli and salmonella bacteria. Brush heads kept in the UV sterilizer after use stays sanitary.



# Aquapick Toothbrush Charcoal & Silver



## Charcoal Toothbrush

**Aquapick Toothbrush Charcoal** is an advanced toothbrush containing a charcoal component and having a function of emitting far-infrared rays and removing the odor with a strong adsorption.




4EA



SUPER SOFT  
Containing charcoal ingredient in the bristle



**SUPER SOFT Fine bristles** effectively eliminates microscopic debris and plaque between teeth and gums because one end of the elastic and smooth of the bristle of 0.18mm is processed to a fine cylindrical shape of 0.03mm or less, special functionality.



Silver

## Silver Toothbrush

**Aquapick Toothbrush Silver** is an advanced toothbrush containing silver component of 99.9% purity and having its strong antibacterial effect proven.



12EA



SUPER SOFT  
Containing silver ingredient in the bristle



# Aquapick Toothbrush Sterilizer AQ-201



Product Name: Aquapick Toothbrush Sterilizer  
Model: AQ-201  
Electrical Currency: AC220V / 60Hz  
Power Consumption: 5W  
Dimension: 240x62x76mm  
Cord Length: 1.3m  
Weight: 270g



Heating part is protected by special Kapton tape, which is effective for preventing electric shock and blocking electric wave.

- Semi-permanently usable without changing the lamp  
Easily sterilize tooth brush which is used everyday.
- 24 hours drying & sterilizing germs  
The PTH heater, which uses infrared radiation, sterilize for 24 hours. This keeps tooth brush clean at all times.

# Aquapick Portable Toothbrush Sterilizer AQ-203



Product Name: Aquapick Portable Toothbrush Sterilizer  
Model: Translucency AQ-202  
Transparency AQ-203  
Longevity of lamp: approximately 20,000 hours  
Dimension: 210x50x25mm  
Weight: 100g (excluding the battery)



CCGL UV Lamp



Usable for 4~6 month



Wall holder



Convenient cap



ventilation



# Aquapick Dental Mask

## AQ-50



Product Name: Aquapick Dental Mask  
Model: AQ-50  
Type: For operation use, For general use  
Amount: 50EA (1BOX)  
Function: Diagnosis, infection prevention  
Storage: store at room temperature  
(1~30℃)  
Dimension: 175x95mm  
Registration Number: 5008

BEF


Effective for blocking **over 99.9%** of bacteria

PEF


Effective for blocking **over 99.9%** of particle

FDA

**100% non-woven fabric** which is harmless



- **FDA** certified non-asbestos material with excellent durability  
| 100% Polypropylene of safe material and passed fire retardant test



- Hypoallergenic, prevents trouble on skin  
| Reduces risk of causing skin trouble because it **doesn't contain any latex, fiberglass, dye, or fluorescent material**



- Fresh air-permeability, soft elasticity  
| **Double wire** provides comfortable wearing sensation  
Outstanding air permeability enables comfortable breathing



- Triple filter, excellent for blood permeate resistance  
| **Blocks bacteria & dust, prevent respiratory disease.**  
Filtering efficiency was tested through laboratories in the US



\* Efficiency certified by the US Nelson Laboratories

# Aquapick Whitening Toothfoam



- Certified product **FDA** & **MFDS**
- Brushing & gargling at once!
- Whitening & Strengthening teeth



It can be used with  
Aquapick oral irrigator.

Fill up the Aquapick water tank with water and pump Aquapick Whitening Toothfoam 2~3 times inside. Use it after diluting the foam with water.

Product name: Aquapick Whitening Tooth foam    Volume: 150ml

Smooth bubbling foam cleaning the teeth!  
Innovative Whitening Toothfoam

- | Fine bubbles
- | Xyllitol for prevention of dental cavities

Hydroxyapatite for strong and white teeth  
(Hydroxyapatite: a major component and an essential ingredient of bone and teeth. It gives bones and teeth their rigidity)

Green tea extract to freshen up your mouth

Good for swollen and bleeding gums  
(Allatoin: prevention of gum disease)

| Contains Lemon oil for shiny teeth

Abrasive-free  
Surfactant-free  
It is good for toothache since it is  
Abrasive-free



- |                  |                       |                |                                  |
|------------------|-----------------------|----------------|----------------------------------|
| Whitening effect | Removal of bad breath | Prevent cavity | Prevent inflammation of the gums |
|------------------|-----------------------|----------------|----------------------------------|



# Aquapick Cleaning Tablet

## AQ-10



Same as actual size

- With one Aquapick cleaning tablet, it takes care of the internal cleaning for oral irrigator!



▶ Oral irrigator for home use



▶ Portable oral irrigator




Product Name: Aquapick Cleaning Tablet  
Model: AQ-10  
Use: For cleaning oral irrigator  
Recommended frequency of use: More than 2~4 times per month  
Storage: Airtight container, room temperature (0~30℃)  
How to use: Put 1 tablet in 200~600ml of water and rinse 1~2 times  
Amount: 12ea/package

### ● Cleaning even invisible spots!



✓ Fur in water tank



✓ Water streaming spots



✓ Substance within the hose



✓ Invisible spots



### ● Great sterilizing effect

Sterilize 99.9% of Escherichia coli, Staphylococcus and Mutans  
- Certificate authority: Korea Testing & Research Institute (KTR)



### ● Made safe

Alcohol-free, Chlorine-free, Dye-free  
Individual packaging facilitates easy use



### ● Sanitary cleaning

Even after cleaning once, it gets as clean as drinking water  
- Certificate authority: Korea Environment & Water Works Institute (KEWWI)



### ● Safe Type 2 Cleaner

Type 2 cleaner is used for kitchenware cleaning  
- Certificate authority: Korea Conformity Laboratories (KCL)

NEW

# Aquapick Electric Kettle

## AQ-400



Product Name: Aquapick Electric Kettle  
Model: AQ-400  
Color: Black, White  
Components: Main body, Lid case,  
Power base.  
Capacity: 0.6L  
Electrical Currency: 220~240V / 60Hz  
Power Consumption : 1100~1310W  
Dimension: 268x134x181mm (WxDxH)  
Weight: 920g



### • Best nozzle

Nozzles designed to pour water neatly and accurately wherever you want creates an optimal stream of hot water wherever it needs to be.



### • One touch power button

Heating starts by pressing the power button located at the bottom of the handle, and automatically shuts off the power when the water temperature reaches 100°C.



### • 360 rotating power base

360 rotating heating plate makes it convenient to hold the handle in any direction, and the bottom of the power base is wrapped around the cord storage and adjustment to fit the space used.







*“Improving Oral Health Globally,”*

Aquapick has received international **approval for safety and quality**



EU certification



Korean certification  
for safe electric  
appliances



Japanese  
certification



European certification for  
meeting the approved  
electrowave standard



Recommended by  
the Korean academy  
of oral health



Korean intellectual  
Property office



FDA certification